

MEXSOL

CITRASOL

BICSOL

LACSOL

Product Portfolio
2024

MEXSOL

Calcium-free dialysis solutions intended for citrate renal replacement therapies

Name	Reference nr.	Na ⁺ mmol/L	K ⁺ mmol/L	Ca ²⁺ mmol/L	Mg ²⁺ mmol/L	Cl ⁻ mmol/L	HPO ₄ ²⁻ mmol/L	HCO ₃ ⁻ mmol/L	Glucose mmol/L
MEXSOL K2 Bi20 P+	722054	133	2	0	1	115.7	1.25	20	7.8
MEXSOL K2 Bi20 P0	722052	133	2	0	1	116.9	0	20	7.8
MEXSOL K4 Bi20 P+	742054	133	4	0	1	117.7	1.25	20	7.8
MEXSOL K4 Bi20 P0	742052	133	4	0	1	118.9	0	20	7.8

Packaging: 5L double-chamber PP bag, 2 bags/box, 96 or 120 bags/pallet

Connectors: Luer Lock + Safe Lock, Nip Port + Luer Lock

Multilanguage IFUs and labelling available

MEXSOL G

Calcium-free dialysis solutions intended for citrate renal replacement therapies

Name	Reference nr.	Na ⁺ mmol/L	K ⁺ mmol/L	Ca ²⁺ mmol/L	Mg ²⁺ mmol/L	Cl ⁻ mmol/L	HPO ₄ ²⁻ mmol/L	HCO ₃ ⁻ mmol/L	Lactate mmol/L	Glucose mmol/L
MEXSOL G K2 Bi22	722250	140	2	0	0.75	120.7	0.8	22	3	6.1
MEXSOL G K4 Bi22	742250	140	4	0	0.75	122.7	0.8	22	3	6.1

Packaging: 5L double-chamber PP bag, 2 bags/box, 96 or 120 bags/pallet

Connectors: Luer Lock + Tulip (Injection Port)

Multilanguage IFUs and labelling available

MEXSOL H

Calcium-containing dialysis solutions intended for renal replacement therapies

Name	Reference nr.	Na ⁺ mmol/L	K ⁺ mmol/L	Ca ²⁺ mmol/L	Mg ²⁺ mmol/L	Cl ⁻ mmol/L	Lactate mmol/L	HCO ₃ ⁻ mmol/L	Glucose mmol/L
MEXSOL H K0 Bi35	703560	140	0	1.5	0.5	109	0	35	5.55
MEXSOL H K0 Bi35 Mg+	703561	138	0	1.5	0.75	108	0	35	5.55
MEXSOL H K1 Bi35	713560	140	1	1.5	0.5	110	0	35	5.55
MEXSOL H K1,5 Bi35	763560	140	1.5	1.5	0.5	111	0	35	5.55
MEXSOL H K2 Bi35	723560	140	2	1.5	0.5	111	0	35	5.55
MEXSOL H K2 Bi35 Mg+	723561	138	2	1.5	0.75	110	0	35	5.55
MEXSOL H K4 Bi35	743560	140	4	1.5	0.5	113	0	35	5.55
MEXSOL H K4 Bi35 Mg+	743561	138	4	1.5	0.75	112	0	35	5.55
MEXSOL H K1 L40	714060	140	1	1.5	0.5	105	40	0	5.55
MEXSOL H K1 L45	714560	140	1	1.5	0.5	100	45	0	5.55
MEXSOL H K1,5 L40	764060	140	1.5	1.5	0.5	106	40	0	5.55
MEXSOL H K1,5 L45	764560	140	1.5	1.5	0.5	101	45	0	5.55
MEXSOL H K2 L40	724060	140	2	1.5	0.5	106	40	0	5.55
MEXSOL H K2 L40 Mg+	724061	140	2	1.5	0.75	107	40	0	5.55
MEXSOL H K2 L45	724560	140	2	1.5	0.5	101	45	0	5.55

Packaging: 5L double-chamber PP bag, 2 bags/box, 96 or 120 bags/pallet, 5L PP bag (lactate solution)

Connectors: Luer Lock + Tulip (Injection Port), Luer Lock + Luer Lock

Multilanguage IFUs and labelling available

BICSOL

Calcium-containing dialysis solutions intended for renal replacement therapies

Name	Reference nr.	Na ⁺ mmol/L	K ⁺ mmol/L	Ca ²⁺ mmol/L	Mg ²⁺ mmol/L	Cl ⁻ mmol/L	HCO ₃ ⁻ mmol/L	Glucose mmol/L
BICSOL 35-108 K0	703550	138	0	1.5	0.75	108	35	5.55
BICSOL 35-110	723550	138	2	1.5	0.75	110	35	5.55
BICSOL 35-112 K4	743550	138	4	1.5	0.75	112	35	5.55

Packaging: 5L double-chamber PP bag, 2 bags/box, 96 or 120 bags/pallet

Connectors: Luer Lock + Injection Port

Multilanguage IFUs and labelling available

LACSOL

Calcium-containing dialysis solutions intended for renal replacement therapies

Name	Reference nr.	Na ⁺ mmol/L	K ⁺ mmol/L	Ca ²⁺ mmol/L	Mg ²⁺ mmol/L	Cl ⁻ mmol/L	Lactate mmol/L	Glucose mmol/L
LACSOL 40-107	541131	140	2	1.5	0.75	107	40	5.55
LACSOL 45-100	541133	140	1	1.5	0.50	100	45	5.55
LACSOL 45-101	541134	140	2	1.5	0.50	101	45	5.55

Packaging: 5L PP bag, 2 bags/box, 96 or 120 bags/pallet

Connectors: Luer Lock + Luer Lock

Multilanguage IFUs and labelling available

CITRASOL 0.5%

Anticoagulation solution intended for citrate renal replacement therapies

Name	PP bag volume	Reference nr.	Na ⁺ mmol/L	Citrate mmol/L	Cl ⁻ mmol/L	Bags/box	Bags/pallet
CITRASOL 0.5%	5000 mL	603050	140	18	86	2	96

Connectors: Luer Lock + Tulip (Injection Port)

Multilanguage IFUs and labelling available

CITRASOL 4%

Anticoagulation solutions intended for citrate renal replacement therapies

Name	PP bag volume	Reference nr.	Na ⁺ mmol/L	Citrate mmol/L	Bags/box	Bags/pallet
CITRASOL 4%	250 mL	601030*	404.6	136.4	18	1512
CITRASOL 4%	250 mL	601031**	404.6	136.4	18	1512
CITRASOL 4%	1000 mL	601010	404.6	136.4	10	480
CITRASOL 4%	2000 mL	601023	404.6	136.4	5	240

Connectors: Luer Lock + Safe Lock

* 601030 connectors: Luer Lock + Luer Lock

** 601031 connectors: Tulip + Injection Port

Multilanguage IFUs and labelling available

CITRASOL ACD-A

Acid Citrate Dextrose – Formula A anticoagulant solution for apheresis methods

Name	PP bag volume	Reference nr.	Na ⁺ mmol/L	Citrate mmol/L	Glucose mmol/L	Bags/box	Bags /pallet	Connector
CITRASOL ACD-A	250 mL	604030	224.4	112.9	123.6	18	1512	Luer Lock
CITRASOL ACD-A	250 mL	604031	224.4	112.9	123.6	18	1512	Infusion Port
CITRASOL ACD-A	500 mL	604050	224.4	112.9	123.6	10	660	Luer Lock
CITRASOL ACD-A	500 mL	604051	224.4	112.9	123.6	10	660	Infusion Port
CITRASOL ACD-A	750 mL	604080	224.4	112.9	123.6	14	672	Luer Lock
CITRASOL ACD-A	750 mL	604081	224.4	112.9	123.6	14	672	Infusion Port

Multilanguage IFUs and labelling available

CITRASOL 4%

4% sodium citrate anticoagulant solution for apheresis methods

Name	PP bag volume	Reference nr.	Na ⁺ mmol/L	Citrate mmol/L	Bags/box	Bags/pallet
CITRASOL 4%	250 mL	601030*	404.6	136.4	18	1512
CITRASOL 4%	250 mL	601031**	404.6	136.4	18	1512

Multilanguage IFUs and labelling available



Manufacturer:

MEDITES PHARMA, spol. s r.o.

1. máje 2635

756 61 Rožnov pod Radhoštěm

Czech Republic

MEDITES PHARMA
spol. s r.o.
Rožnov pod Radhoštěm

www.meditespharma.cz